



Single-cell RTq-PCR: Raw data treatment

Arnaud Bonnaffoux

Single-Cell day




Équipe DRACULA




association nationale
recherche technologie

INTRODUCTION



(Richard et al. (2016) PLoS Biology)





SOURCE OF TECHNICAL NOISE



Raw data treatment/checks steps

1. Deletion of abnormal amplifications (check qPCR)
2. « Black cells » deletion (check drop-outs)
3. Deletion of low-detected gene (check RT)
4. Spike normalisation (RT bias correction)
5. Conversion of Ct in molecule counts
6. Correction factor (dilution bias correction)
7. Gene triplicate (check qPCR)



1- Deletion of abnormal amplified cell/gene



- Export Ct data in .csv file
- Observe « manually » melt curve of one gene for all cells using Fluidigm software
- If a pic is shifted ($\Delta > 0.6^\circ$) compared to main pic, set Ct value for couple cell/gene to UD (UnDefined)

2- « Black » cells deletion


- « Black » cells = Spikes with low Ct, all other genes have important Ct ($Ct > 30 \neq \text{UnDefined}$)
- Even if genes expression is bursty (important probability of null value), the probability to have ALL genes to 0 is very low
- Proportion of black cells can be important :
 - AR11 : 27%
 - AR78 : 54%
 - AR85 : 32%
- Check absence of « black » cells using sum of Ct per cell



AR11

3- Deletion of low-detected gene

- For genes with important UnDefined values ($\neq 0$) in cells ($> 80\%$), sample size is too low to be considered
- Only few genes are concerned
 - AR78 : 3 / 93
 - AR85 : 2 / 93



4- SPIKE : check RT + normalisation

- At least one spike shall be detected and Ct < 30, else cell is deleted
- If spikes correlation is OK (R>0,7) normalise Ct values of a cell with its spike values
- Spike correlation is « stochastic » (???)

$$C_{q_{i,j}} = \widehat{C_{q_{i,j}}} - (\overline{C_{q_i}} - \overline{C_{q_0}})$$


Raw Ct

All cells mean spikes Ct

Ct normalized


Spike Ct

spike correlation with Cor_coeff = 0.91




AR11

spike correlation with Cor_coeff = 0.25



AR78

spike correlation with Cor_coeff = 0.77



AR85


5- Conversion of Ct in molecule count

$$m_{i,j} = 96 \times 45 \times 2^{30-22-Cq_{i,j}}$$


- 30 = detection threshold
- 22 = cycles of pre-amplification
- 96 = dilution in 96 wells
- 45 = sampling factor for qPCR measurement

6- Correction factor: Total distribution centering


- To correct possible bias due to imprecision of dilution factor or detection threshold compare your total gene expression distribution (in all cells) with literature distribution




Schwanhausser et al., 2013, Nature



AR78 (*1)






AR85 (*1)



AR85 (*2)

7-GENE TRIPLICATE : check qPCR

- Use 3 identical primers for one gene (RPL22L1)
- Usually qPCR replication is quite perfect



Raw data treatment/checks steps

1. Deletion of abnormal amplifications (check qPCR)
2. « Black cells » deletion (check drop-outs)
3. Deletion of low-detected gene (check RT)
4. Spike normalisation (RT bias correction)
5. Conversion of Ct in molecule counts
6. Correction factor (dilution bias correction)
7. Gene triplicate (check qPCR)



Olivier Gandrillon
Angélique Richard
Elodie Vallin
Sandrine Giraud
Ulysse Herbach



DRACULA group